

Data Acquisition Tutorial Qcodes

April 2, 2024

```
[1]: %matplotlib inline
import numpy as np
from pathlib import Path
import qcodes as qc
import matplotlib.pyplot as plt

## Multidimensional scanning module
from qcodes.dataset import (
    Measurement,
    initialise_or_create_database_at,
    load_by_guid,
    load_by_run_spec,
    load_or_create_experiment,
    plot_dataset,
)

## Dummy instruments for generating synthetic data
from qcodes.instrument_drivers.mock_instruments import (
    DummyInstrument,
    DummyInstrumentWithMeasurement,
)

## Using interactive widget
from qcodes.interactive_widget import experiments_widget
```

0.1 Qcodes Data Accumulation Tutorial

The below commands are for setting up “dummy instruments.” This means that they are simulators of instruments, in this case a signal generator and a multimeter, that have preprogrammed behaviors. This document is designed to use these dummy instruments to simply demonstrate the data measurement architecture and data saving we will use in QcOdEs. You don’t need to fully understand these dummy instruments themselves, just the measurement structures.

```
[2]: # A dummy signal generator with two parameters ch1 and ch2
dac = DummyInstrument('dac', gates=['ch1', 'ch2'])

# A dummy digital multimeter that generates a synthetic data depending
# on the values set on the setter_instr, in this case the dummy dac
dmm = DummyInstrumentWithMeasurement('dmm', setter_instr=dac)
```

```
dac.ch1(1.1)
dmm.v1()
```

[2]: 3.9097352904646567

The cell below is used to initialize an experiment for saving and reaccessing experiment data as well as assigning parameters for data saving.

```
[3]: #We start by stopping the laser pulsing. This way we can properly initialize.
initialise_or_create_database_at(Path.cwd() / "Experiment Name.db") #Experiment_
↳name for data saving file to open or create
experiment = load_or_create_experiment(
    experiment_name='Experiment Name',
    sample_name="" #Can use for sample names, but we don't vary the sample, so_
↳can be left blank
)
#Use this to make a parameter out of anything!
DAC = qc.ManualParameter('Setting', unit='V')
DMM = qc.ManualParameter('Signal', unit='V')
meas = Measurement(exp=experiment, name='Experiment Name')
meas.register_parameter(DAC) #register the first independent parameter
meas.register_parameter(DMM)
```

There is (are) already experiment(s) with the name of Experiment Name and sample name of some_sample in the database.

[3]: <qcodes.dataset.measurements.Measurement at 0x1e50832e200>

```
[4]: dac_range = np.linspace(0, 10, 101)
with meas.run() as your_data: #This begins the experiment and assigns the name_
↳of your dataset.
    for dac_val in dac_range: #Your usual for loop to get a voltage to set the_
↳DAC to.
        #####
        #Here you put any actions needed for each step of the experiment
        dac.ch1(dac_val)
        #####
        #This is used for adding your results to your data set in the form_
↳(parameter, value).
        your_data.add_result((DAC, dac_val),
                             (DMM, dmm.v1()))
        output_data = your_data.dataset #Though you can directly use the variable_
↳your_data, if running several experiments, just note you can access the_
↳dataset values
                                     #By using this command.
```


Starting experimental run with id: 10.

Once this is finished, you can convert this to a familiar format and do your usual data analysis as shown below. Note that, depending on how you

```
[5]: data = output_data.to_pandas_dataframe()
      print(data)
      plt.plot(data["Setting"], data["Signal"])
      plt.xlabel('DAC Setting(V)')
      plt.ylabel('DMM Signal(V)')
      plt.title('Ch1 DAC setting vs DMM signal graph')
      plt.show()
```

	Setting	Signal
0	0.0	5.026318
1	0.1	4.869901
2	0.2	4.815182
3	0.3	4.728395
4	0.4	4.743018
..
96	9.6	0.824237
97	9.7	0.752472
98	9.8	0.730891
99	9.9	0.563374
100	10.0	0.648558


[101 rows x 2 columns]



0.2 Retrieving Data

Data retrieval is very simple in this software. You are given an experimental run ID with every new run in a given experiment. To access the data, you can run the code block below with the desired run ID.

```
[6]: loaded_data = load_by_run_spec(captured_run_id=9) #Substitute 1 for the desired_
      ↪experiment number
      panda_loaded_data = loaded_data.to_pandas_dataframe() #Converts data into a_
      ↪panda data frame.
      plt.plot(panda_loaded_data["Setting"], panda_loaded_data["Signal"])
      plt.xlabel('DAC Setting(V)')
      plt.ylabel('DMM Signal(V)')
      plt.title('Ch1 DAC setting vs DMM signal graph')
      plt.show()
```



0.3 Multiple Variables and Contour Maps

You can also do measurements with multiple variables.

```
[7]: #We start by stopping the laser pulsing. This way we can properly initialize.
initialise_or_create_database_at(Path.cwd() / "Experiment Name 2.db")
    ↪ #Experiment name for data saving file to open or create
experiment = load_or_create_experiment(
    experiment_name='Experiment Name 2',
    sample_name="" #Can use for sample names, but we don't vary the sample, so
    ↪ can be left blank
)
#Use this to make a parameter out of anything!
DAC1 = qc.ManualParameter('Setting1', unit='V')
DAC2 = qc.ManualParameter('Setting2', unit='V')
DMM = qc.ManualParameter('Signal', unit='V')
meas = Measurement(exp=experiment, name='Experiment Name 2')
meas.register_parameter(DAC1)
meas.register_parameter(DAC2)
meas.register_parameter(DMM)
```

There is (are) already experiment(s) with the name of Experiment Name 2 and sample name of some_sample in the database.

[7]: <qcodes.dataset.measurements.Measurement at 0x1e50a2b7220>

```
[8]: dac_range1 = np.linspace(-1, 1, 101)
dac_range2 = np.linspace(-1, 1, 101)
with meas.run() as your_data: #This begins the experiment and assigns the name
    of your dataset.
    for dac1_val in dac_range1:

        for dac2_val in dac_range2: #Your usual for loop to get a voltage to
            set the DAC to.

                #####
                #Here you put any actions needed for each step of the experiment
                dac.ch1(dac_val)
                dac.ch2(dac_val)
                #####
                #This is used for adding your results to your data set in the form
                (parameter, value).
                your_data.add_result((DAC1, dac1_val),
                                    (DAC2, dac2_val),
                                    (DMM, dmm.v1()))


        output_data = your_data.dataset

graph_data = output_data.to_pandas_dataframe()
fig=plt.figure()
dac_ch1_dim=np.linspace(dac_range1[0], dac_range1[-1], len(dac_range1))
dac_ch2_dim=np.linspace(dac_range2[0], dac_range2[-1], len(dac_range2))
dac1, dac2 = np.meshgrid(dac_ch1_dim, dac_ch2_dim)
sig = np.add(list(graph_data["Signal"]), 0)
sig = sig.reshape(len(dac_ch1_dim),len(dac_ch2_dim)).T
plt.contourf(dac1, dac2, sig, 500) #The number chosen will give differeing

fig.suptitle('Ch1, Ch2 DAC for DMM signal Contour Map of Signal')
plt.xlabel('Ch1 (V)')
plt.ylabel('Ch2 (V)')
plt.colorbar()
#plt.pcolormesh(dac1, dac2, sig, shading='auto', cmap='Blues') #Can be used to
change the color
fig.canvas.draw_idle()
```

Starting experimental run with id: 40.

Ch1, Ch2 DAC for DMM signal Contour Map of Signal



The data here is simply just noise, so there isn't any point reading into it. Following a procedure similar to this, you can acquire data in such a way that you can get 3D graphs or color graphs similar to the one above.